


JorVet™

Gastric Lavage Kit # J-115


A complete kit designed for emergency treatment of gastric dilatation/volvulus (Bloat) and acute poisonings. This include 3 different sizes of large bore, semi-rigid, lavage tubes that pass down the esophagus and into the stomach. The large size tube allows the passage of a stomach irrigation system that will flush out ingesta from the upper GI. The metal marker in the tip of the stomach tube can be checked for location by radiograph.



J-115

Kit Includes:

- 1) Semi-rigid polyethylene tube, 1 1/2" OD; 1 1/4" ID; 48" long
- 2) Semi-rigid polyethylene tube, 1" OD; 7/8" ID; 42" long
- 3) Semi-rigid polyethylene tube, 3/4" OD; 5/8" ID; 36" long
- 4) Dura Clamp (to prevent backflow)
- 5) Reducing Clamp
- 6) Faucet adapter
- 7) PVC tubing, 1/2" OD; 3/8" ID x 100cm
- 8) 23 x 76cm stomach tube with metal x-ray marker



Jorgensen Laboratories, Inc.

1450 Van Buren Ave.

Loveland, CO 80538

(800) 525-5614

www.JorVet.com • info@JorVet.com