



**NEW!**  
**Piezo Electric Scaler Polisher Combo Unit**  
**from JorVet**

We now offer a combo unit with the very popular Piezo hand piece with a sleek design at an affordable price!



- Contemporary design with a small footprint
- Operates at 29kHz
- Includes in-line water filter, 3 scaler tips (2 spatula, 1 universal), micro-motor and polisher hand-piece, prophy angle, prophy cups and paste starter kit.
- Comes with quick release attachment for water tank set-up (optional water tank J-452d4)
- Micromotor variable speed 2,000-30,000rpm, with forward/reverse
- Large, easy to read LED power display
- Foot pedal operation
- Optional drilling and sectioning attachments available
- **6 year** limited warranty (90 days on tips and prophy angle, 2 years on handpiece)
- Dimensions: 6"L x 8"W x 3.5"H
- Weight 6.5 lbs



**J-841x**  
**Piezo Electric Scaler Polisher Combo Unit**

**Jorgensen Laboratories, Inc.**

1450 Van Buren Avenue  
Loveland, CO 80538  
(800) 525-5614

[www.JorVet.com](http://www.JorVet.com) [Info@JorVet.com](mailto:Info@JorVet.com)